
ClariFix® Cryotherapy device

from stryker

FOR CHRONIC RHINITIS – RUNNY NOSE, CONGESTION, OR POST-NASAL DRIP



Are you tired of always having a runny nose, nasal congestion, or postnasal drip? For the 87 million Americans who suffer from <u>chronic</u> <u>rhinitis</u>. Stryker's ClariFix device goes directly to the root of chronic rhinitis– the overactive nerves. Using a minimally invasive treatment called cryotherapy, the ClariFix device applies controlled freezing to target tissue under direct visualization.



HOW CLARIFIX WORKS

ClariFix cryotherapy works by creating second degree

nerve damage which has demonstrated rapid regeneration of mucosal structures and function. Additionally, large blood vessels have shown to be resistant to the freezing and their function is not impaired.

CLARIFIX CRYOTHERAPY HAS BEEN SHOWN TO BE SAFE AND EFFECTIVE.

- BACKED BY 5 CLINICAL STUDIES
- MINIMALLY INVASIVE
- SIDE EFFECTS MAY BE MILD, FLEETING, AND RESOLVE ON THEIR OWN
- More than 50,000 patients treated since 2017
- MINIMAL DOWNTIME AND RECOVERY
- CAN BE PERFORMED UNDER LOCAL ANESTHESIA